medartis®

M-0310_v0 - 05/2010

Medartis®

Medartis®

Instructions for Use:	MODUS TTS Tooth Splint, Titanium (Titanium Trauma Splint)
Gebrauchsanweisung:	MODUS TTS-Zahnschiene, Titan (Titanium Trauma Splint)
Mode d'emploi :	Attelle MODUS TTS, Titane (Titanium Trauma Splint)
Istruzioni d'uso:	MODUS TTS ferula dentale, titanio (Titanium Trauma Splint)

Manufacturer/Hersteller/Fabricant/Produttore:
Medartis AG, Hochbergerstrasse 60E, CH-4057 Basel (Switzerland/Schweiz/Suisse/Svizzera)
Phone/Telefon/Téléphone/Telefono +41 61 633 34 34
Telefax/Telefax/Télécopie/Telefax +41 61 633 34 00
E-mail info@medartis.com

Distribution/Vertrieb/Distribution/Distributore:
Please contact your national distributor or Medartis AG, Basel.
Bitte kontaktieren Sie Ihre lokale Vertretung oder die Medartis AG, Basel.
SVP veuillez contacter votre distributeur local ou Medartis AG, Bâle.
Per favore mettersi in contatto con il vostro distributore nazionale o Medartis AG, Basilea.

English
Introduction
These instructions for use are for a product line of Medartis AG, 4057 Basel/Switzerland. All of the directions provided in these instructions must be followed.

Notes Regarding the Delivered Items
The TTS Tooth Splint may only be accepted when the manufacturer’s labels and packaging are undamaged and unopened at the time of delivery. If this is not the case, the questionable items must be returned to Medartis AG, Basel, Switzerland or to the relevant Medartis branch or distribution partner within ten working days.
The TTS Tooth Splint is only intended for a single use and is not designed to be reused. It is delivered in a **NON-STERILE** condition and must be appropriately prepared before initial use. All packaging must be removed before preparation.

Product Material
The TTS Tooth Splint consists of pure titanium (ASTM F67, ISO 5832-2). The titanium material used is biocompatible, corrosion-resistant and non-toxic in a biological environment.

Art. No.	Description	Size
TTS- 1002	TTS Tooth Splint	52 x 0.2 mm (length x thickness)
TTS- 1004	TTS Tooth Splint	100 x 0.2 mm (length x thickness)

Proper Use
The TTS Tooth Splint is used for the treatment of traumatized or replanted teeth.

Indications
TTS Tooth Splint is used for:

- Stabilization after repositioning or replantation (acute cases)
- Stabilization after intentional (planned) replantation

Contraindications
None are known.

Possible Complications

- Irritation (inflammation) of the inner part of the lip. The TTS Tooth Splint must be fixed so that it lies flat, directly on the surface of the enamel.
- Irritation or inflammation of the gum, caused by placing the TTS Tooth Splint too close to it.
- Speech impairment.
- Use of excessive amounts of composite may cause hygiene problems.
- If adhesion of the composite to the surface of the enamel is inadequate (insufficiently dry, no bonding) the TTS Tooth Splint can become loose. Consequence: **danger of swallowing or aspiration of the splint.**

Warnings and Precautionary Measures

- The product may only be used by medical personnel who hold relevant qualifications.
- As the manufacturer, Medartis recommends that the user read all of the available documents before first use and contact other users who have practical experience with this type of treatment.
- Never use products that have been damaged by transport, improper handling in the hospital, or in any other way!
- The TTS Tooth Splint is intended for single use and may not be reused under any circumstances.
- When used in the mouth, the product must be secured against aspiration.
- Observe the necessary care when using and storing the products:
 - Damage to (such as from improper cutting or bending) and/or scratches on the product can substantially impair the strength of the product and lead to premature breakage.
- The TTS Tooth Splint has been developed and manufactured for a specific purpose and may not be modified by the user.

Instructions for Selecting the Appropriate Product
As the manufacturer, Medartis does not recommend specific surgical procedures for specific patients. The dentist is responsible for choosing the appropriate product for the specific case. Assessment of the patient and correct selection and positioning are prerequisites for clinical success.
The treating dentist should become thoroughly familiarized with the procedure, for example by:

- Carefully studying all the product documentation
- Carefully reviewing the current professional literature
- Consulting with colleagues experienced in this field and with the use of this tooth splint
- Practice handling the tooth splint

Clinical Use of the TTS Tooth Splint
Inserting the Tooth Splint

- Take the TTS Tooth Splint out of its packaging
- Determine the splint length on the patient
- Cut the splint to the required length (e.g. by using a crown-cutter or side-cutter)
- Adapt the splint to the curvature of the teeth (by hand, no instruments required)
- Clean the enamel surfaces by rinsing and drying them
- Etch the enamel surfaces (only at those points on the surfaces where adhesion is intended, as defined by the splint design)
- Apply the adhesive agent (bonding) and carry out light polymerization
- Position the splint by hand or using tweezers
- Fill the hollow, diamond-shaped spaces in the splint to be fixed to the tooth surfaces, using free-flowing composite, then carry out light polymerization
- Check for possible irritation of lips and soft tissues of the cheek by the splint; if necessary modify the splint (re-cutting), especially at its ends
- Splinting time, max. 2-4 weeks

Removing the Splint

- Remove the composite down to the splint level, using rotary instruments (**not up to the enamel surface!**)
- Take hold of the splint at one end and pull it off the crown surfaces like a piece of film
- Remove any composite which remains, taking care not to harm the enamel

Duration of Application for the TTS Tooth Splint
Splinting time: max. 2-4 weeks

Documentation and Traceability
Medartis recommends full clinical and radiological documentation. Each TTS Tooth Splint can be traced back with the aid of its Item No. and Lot No. The product label includes all relevant data for notation on the patient’s record.

Additional Information
Additional information on the product can be requested from your local Medartis branch or distributor. In addition, all the relevant information can be found on the internet at www.medartis.com.

Instructions Regarding Cleaning, Disinfection and Sterilization
The TTS Tooth Splint is delivered **NON-STERILE** and must be cleaned, disinfected and sterilized before each use. This also holds true for initial use after the delivery (cleaning and disinfection after removing the protective transport packaging, and sterilization after being placed in the sterilization container). TTS Tooth Splints which come in contact with patient’s blood or other bodily fluids must be discarded. Thorough cleaning and disinfection are essential for effective sterilization. It is your responsibility to ensure that the TTS Tooth Splint is completely sterile when used, to employ appropriately validated cleaning, disinfecting and sterilizing procedures, to regularly service and inspect the employed devices (disinfector and sterilizer), and to ensure that the validated parameters are maintained for each cycle. The statutory regulations applicable for your country and the hospital’s hygiene requirements must also be observed. This applies in particular to the various instructions for effectively deactivating prions.

Basic Instructions
If possible, use an automated method (disinfector) for cleaning and disinfection. Do not use a manual procedure even with an ultrasonic bath due to the significantly reduced efficiency and potential damage.

Manual Cleaning and Disinfection
Observe the following when choosing cleansers and disinfectants:

- That they are suitable for cleaning and disinfection
- That the cleanser (if applicable) is suitable for ultrasonic cleaning (no foam formation)
- Use a disinfectant with proven effectiveness (such as approval by VAH/DGHH or the FDA, or a CE mark)
- That the employed chemicals are compatible

Manual Cleaning Process

- Place the product in the cleaning bath (the product must be completely covered, the specified reaction time must be observed, and the individual product may not touch an-other product)
- Clean using a soft brush
- Rinse at least twice thoroughly under flowing water (demineralized, at least tap water quality in terms of the bacterial level)
- Inspect the product (see the section "Inspection")

Manual Disinfection

- Place the cleaned and inspected product in the disinfection bath (the product must be completely covered, the specified reaction time must be observed, and the product may not touch another product)
- Rinse at least twice thoroughly under flowing water (demineralized, at least tap water quality in terms of the bacterial level)
- The product must be completely dried immediately afterwards (it is recommendable to dry them using compressed air)
- Inspect the product (see the section, "Inspection")
- Pack the product immediately after it has fully dried

Automated Cleaning and Disinfection
When selecting the disinfector, make sure:

- That the effectiveness of the disinfector is proven (such as approval by VAH/DGHH or the FDA, or a CE mark).
- That the cleaning process includes the following phases in accordance with EN ISO 15883:

Phase	Temperature	Duration	Action
Cleaning	Heat to 93°C (199.4°F)	Device-specific	The cleaner is dispersed and suspended
Thermal disinfection (Ao value >3000)	93°C (199.4°F)	10 min.	Do not add additional cleanser
Rinsing	-	Device-specific	Rinse with demineralized water

Inspection of the TTS Tooth Splint
Before placing the TTS Tooth Splint in the sterilizing container, check it for damages to the surface and contaminants. Separate damaged or non clean TTS Tooth Splints.

Sterilization
Only use the following sterilization methods. Other sterilization methods are not allowed.

Steam Sterilization
All **NON-STERILE** products can be sterilized in an autoclave (EN 13060 and EN 285). For both initial and subsequent sterilization, the following parameters were verified by the manufacturer in accordance with the requirements of the current sterilization standard:

Procedure	Fractionated vacuum method	Flow method
Duration of exposure	≥ 5 min.	≥ 15 min.
Temperature	134°C (273°F)	134°C (273°F)
Drying time	> 20 min. – 30 min.	> 15 min. – 30 min.

In principle, Medartis recommends using the fractionated vacuum method for sterilization with an exposure time of ≥ 18 minutes. Steam sterilization with the gravitation method must be verified with additional validation for the specific product, sterilizer and procedure.In addition, do not use hot-air sterilization, irradiation sterilization, sterilization with formaldehyde or ethylene oxide, and do not use a substitute procedure for sterilizing thermolabile goods such as plasma or peroxide sterilization for the TTS Tooth Splint.Flash sterilization is not recommended. If this method is used, observe the respective country-specific laws, standards, guidelines and instructions. The user is responsible for reviewing these requirements and procuring the related information. After sterilization, the product must be stored in a dry environment.

Reusability
The TTS Tooth Splint may only be brought into contact with an individual patient. TTS Tooth Splints which have come into contact with the blood or bodily fluids of a patient may not be reused. No liability is assumed by the manufacturer in case of non-observance.
If TTS Tooth Splints come in contact with pathogens that are difficult to identify such as varia-tions of Creutzfeldt-Jakob’s disease (confirmed or suspected pathogen), they must be thrown away and may not be reused.

Explanation of Symbols
See reverse

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