

CASE REPORT



Fixation of Ankle Fractures with APTUS 2.8 Fibula Plates and APTUS 3.5 Tibia T-Plates

The Surgeon

Surgeon Name: Dr. T. Schepers MD PhD

Hospital: Amsterdam UMC location Meibergdreef

Dr. Schepers is a trauma-surgeon with over ten years of experience in complex foot ankle injury. He has published over 150 scientific papers on the subject.

Introduction

Ankle fractures are amongst the most commonly encountered injuries. In case of instability these fractures need surgical stabilization. Key to success are anatomical reduction and preventing wound complications.

The Case



Patient Profile

A 36 year old female patient arrived at our Emergency Department by ambulance. There was no relevant prior medical history no relevant medical history. She fell from her scooter and noticed the abnormal position of her ankle.


At the emergency department a radiograph was made, showing a Supination External Rotation Type 4 ankle fracture with a dislocation of the ankle joint. The ankle was promptly reduced and a CT scan was made. Unfortunately the CT-scan still showed a sub-optimal reduction. For which a second reduction was performed and radiographs were made.





Clinical Findings/Preoperative analysis

The CT-scan showed a Weber-B type fibula fracture and a posterior malleolar fracture Bartonicek type 2. On the medial side the deltoid ligament was torn.




Surgical treatment

Seven days after trauma she was operated on in a prone position. A posterolateral approach was performed to address the posterior malleolus. The sural nerve was identified and spared throughout the procedure. An anatomical pre-contoured 3.5 mm T-plate was used to maintain the reduction. Via the same incision the fibula was approached and a 2.8 plate was used following reduction of the fracture.




Intraoperative findings

The APTUS 3.5 T-plate fitted well posteriorly. The APTUS 2.8 Fibula plate allowed for an excellent stabilization of the fracture.



Postoperative treatment

Post-operatively the patient was kept in a cast, non-weight-bearing for 4 weeks and weight bearing for another 2 weeks after which the cast was removed and she was allowed to start weight bearing in a well fitted shoe. She was last seen for follow-up after six months and experienced no pain and had an excellent range of motion.





Conclusion:

Open reduction and internal fixation of an unstable ankle fracture using an anatomically fitted posterior 3.5 Tibia T-plate and a 2.8 Fibula plate provides good results.



References

- 1) Dingemans SA, Lodeizen OA, Goslings JC, Schepers T. Reinforced fixation of distal fibula fractures in elderly patients; A meta-analysis of biomechanical studies. Clin Biomech (Bristol). 2016;36:14-20. doi:10.1016/j.clinbiomech.2016.05.006
- 2) Schepers T, Van Lieshout EM, De Vries MR, Van der Elst M. Increased rates of wound complications with locking plates in distal fibular fractures. Injury. 2011;42(10):1125-1129. doi:10.1016/j.injury.2011.01.009

Disclaimer: This case report presents the outcome of an individual patient case and it does not imply any guarantee or warranty in regard to treatment success. A surgeon must always rely on her or his own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. The professional must always comply with the individual product's Instructions For Use (www.medartis.com/documentation/instructions-for-use) as well as all laws and regulations. Medartis is not giving any medical advice. The devices may not be available in all countries due to registration and/or medical practices. All content – such as texts, video and pictures – was created by healthcare professional mentioned in the case report. For further questions, please contact your Medartis representative (www.medartis.com). This information contains products with CE and/or UKCA marking.

For US only: Federal law restricts this device to sale by or on the order of a physician.

© Medartis 2024. Everything herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Medartis or its affiliates unless otherwise indicated. It is forbidden to redistribute, duplicate or disclose, anything herein, in whole or in part, without the prior written consent of Medartis.