


CASE REPORT



Treatment of acquired flatfoot deformity with TMT instability and HV deformity by performing TMT fusion and HV osteotomy

The Surgeon

Mario Herrera-Perrez

Dr. Mario Herrera is a Doctor of Medicine and Surgery from the University of La Laguna and a reference in advanced foot and ankle surgery, as well as arthroscopic surgery. He also has experience in total ankle prosthesis.

He combines his medical work with teaching as associate professor of Traumatology and Orthopaedics at the Faculty of Medicine of the University of La Laguna. He is also a reviewer for the Spanish Journal of Orthopaedic Surgery and Traumatology (RECOT) and a member of the Editorial Committee of the Spanish Journal of Foot Surgery.

He is currently a specialist at the CECOTEN Centre and is head of section of the Functional Foot and Ankle Unit at the University Hospital of the Canary Islands.

The Case



Patient Profile

A 65 year-old woman was diagnosed of an adult acquired flatfoot deformity with tarsometatarsal instability and hallux valgus deformity (Figs. 1 and 2). After conservative treatment with physiotherapy and insoles, a surgical procedure was considered.



Fig. 1: Hallux valgus deformity.



Fig. 2: Disrupted Meary line in the AAFD.

Figure 1

Figure 2



Surgical Treatment

Through a medial approach of the first tarsometatarsal (TMT) joint, we performed an osteotomy of the base of the M1 parallel to the TMT joint and a second cut of the C1 just parallel to the axis of the M2. After this, we placed the TMT Medartis 2.8 Lapidus Medial plate. Before fixing the plate with screws, we put an interfragmentary compression screw from M1 to C1. In a second step we addressed the hallux valgus (HV) deformity: exostectomy, adductor tenotomy and akin osteotomy were performed.



Postoperative treatment

The aftercare treatment comprised a partial weight bearing in a stiff soled shoe for 6 weeks.



Evolution

After one year follow up, all the parameters are corrected, the foot is well aligned and the Lapidus procedure shows a complete TMT fusion (Figs. 3 and 4).



Fig. 3. TMT fusion and HV procedure.



Fig. 4. Lateral X-ray of the same case showing a normal Meary line.



Conclusion

The modified Lapidus procedure is considered a useful tool for correction of multiplanar deformities in the hallux valgus, including first metatarsal pronation. It offers a greater power of correction compared to most other osteotomies. The combination of a dorsal plate and an interfragmentary screw offers high stability in all planes.



References

- 1) Zambelli R, Baumfeld D. Intraoperative and Postoperative Evaluation of Hallux Valgus Correction: What Is Important?. *Foot Ankle Clin.* 2020;25(1):127-139. doi:10.1016/j.fcl.2019.10.007
- 2) Jagadale VS, Thomas RL. A Clinoradiological and Functional Evaluation of Lapidus Surgery for Moderate to Severe Bunion Deformity Shows Excellent Stable Correction and High Long-Term Patient Satisfaction. *Foot Ankle Spec.* 2020;13(6):488-493. doi:10.1177/1938640019890716
- 3) Haas Z, Hamilton G, Sundstrom D, Ford L. Maintenance of correction of first metatarsal closing base wedge osteotomies versus modified Lapidus arthrodesis for moderate to severe hallux valgus deformity. *J Foot Ankle Surg.* 2007;46(5):358-365. doi:10.1053/j.jfas.2007.05.008
- 4) McAleer JP, Dayton P, DeCarbo WT, et al. A Systematic Approach to the Surgical Correction of Combined Hallux Valgus and Metatarsus Adductus Deformities. *J Foot Ankle Surg.* 2021;60(5):1048-1053. doi:10.1053/j.
- 5) Wagner E, Wagner P. Metatarsal Pronation in Hallux Valgus Deformity: A Review. *J Am Acad Orthop Surg Glob Res Rev.* 2020;4(6):e20.00091. Published 2020 Jun 15. doi:10.5435/JAAOSGlobal-D-20-00091

Disclaimer: This case report presents the outcome of an individual patient case and it does not imply any guarantee or warranty in regard to treatment success. A surgeon must always rely on her or his own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. The professional must always comply with the individual product's Instructions For Use (www.medartis.com/documentation/instructions-for-use) as well as all laws and regulations. Medartis is not giving any medical advice. The devices may not be available in all countries due to registration and/or medical practices. All content – such as texts, video and pictures – was created by healthcare professional mentioned in the case report. For further questions, please contact your Medartis representative (www.medartis.com). This information contains products with CE and/or UKCA marking.
For US only: Federal law restricts this device to sale by or on the order of a physician.

© Medartis 2024. Everything herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Medartis or its affiliates unless otherwise indicated. It is forbidden to redistribute, duplicate or disclose, anything herein, in whole or in part, without the prior written consent of Medartis.