



CASE REPORT

Recalcitrant Scaphoid Nonunion Treated with a Palmar Angular-Stable Scaphoid Plate and Iliac Crest Cancellous Bone Grafting

The Surgeon

PD Dr. Stefan Quadlbauer, FEBHS
AUVA Trauma Hospital Lorenz Böhler, Vienna, Austria

Dr. Quadlbauer is a specialist in orthopedic and trauma surgery with sub-specialization and board certification in hand surgery at the Austrian Workers' Compensation Board (AUVA) Trauma Hospital Lorenz Böhler in Vienna. He is the head of the hand surgery team and is a member of the teaching faculty of the Vienna Hand Courses.

His clinical expertise encompasses complex traumatic and reconstructive conditions of the hand and upper extremity. Research activities include clinical studies on surgical outcomes, innovations in hand and trauma surgery, and the optimization of rehabilitation strategies.

Introduction

Scaphoid fractures account for approximately 10% of all hand fractures¹. Management of scaphoid fractures that fail to heal presents a significant challenge for hand surgeons. Particularly, recalcitrant scaphoid nonunions—defined by prior failed surgery, avascular necrosis, bone loss exceeding 7 mm or nonunion duration of more than two years—are very difficult to treat and are associated with a substantially reduced probability of union. In such cases, salvage procedures, including wrist denervation, radial styloidectomy, four-corner arthrodesis or proximal row carpectomy are frequently considered. However, palmar angular-stable scaphoid plating has emerged as a reliable surgical option, providing stable fixation and promoting bone consolidation even in these complex cases¹⁻⁷.

The Case



Patient History / Profile

A 22-year-old male sustained a wrist injury during a football game six years prior to the surgical intervention. Since the initial trauma the patient reported persistent wrist pain. He was subsequently referred to our outpatient hand clinic, where imaging studies confirmed a scaphoid nonunion. The patient had not undergone any prior surgical treatment for this injury.



Preoperative Clinical Findings

The patient exhibited an unrestricted range of motion in the wrist but reported pain during load-bearing activities. Radiographs (Figure 1a–1d) and computed tomography (CT) scans (Figure 1e, 1f) demonstrated a scaphoid nonunion with a severe humpback deformity and dorsal intercalated segment instability (DISI) of the lunate⁸. Advanced cystic changes and sclerosis were present at the nonunion site. Magnetic resonance imaging (MRI) revealed bone marrow edema in both fragments (Figure 1g) and absence of fat-suppression signal in the proximal scaphoid pole (Figure 1h), consistent with avascular necrosis. No significant osteoarthritic changes were observed in the wrist.




Figure 1: Pre-OP imaging of the scaphoid: a-d radiographs; e,f CT images; g,h MRI images



Surgical Treatment

The scaphoid nonunion was approached through a palmar incision over the flexor carpi radialis (FCR) tendon. After preparation of the FCR tendon, the tendon's sheath was incised, and the radioscapocapitate ligament identified and longitudinally split to expose the scaphoid¹¹. A palmar cortical window was created, and fibrous and necrotic tissue was debrided. To correct the humpback deformity and restore scaphoid length, the DISI position of the lunate was reduced using the Linscheid maneuver. In this maneuver, the wrist was flexed under image intensifier control, bringing the lunate into neutral position. Temporary transfixion of the lunate to the radius was achieved using a 1.4 mm K-wire (Figure 2a, 2b). An additional 1.2 mm K-wire was inserted perpendicularly into the distal scaphoid fragment and used as a joystick to facilitate final reduction.

Following restoration of scaphoid length and shape, residual sclerotic and fibrous tissue was removed using a low-speed spiral burr. The reduction was temporarily stabilized with a transfixing K-wire, and both fragments were additionally drilled until punctate bleedings were observed. The defect was filled with cancellous bone graft harvested from the iliac crest and impacted.

Due to the anatomy of the patient and size of the nonunion defect, the large version of the scaphoid plate was selected and positioned under image intensifier control. The plate was prefixed distally with a K-wire inserted through a screw hole. After verifying the correct plate position, locking screws were inserted after drilling and assigning the screw length (Figure 2c–2g).




Figure 2: Intra-OP images: a,b Linscheid maneuver to correct the DISI position of the lunate; c-f radiographs of plate position; g photograph of plate in situ

One session of high-energy extracorporeal shockwave therapy (ESWT) was applied immediately postoperatively, delivering 3,000 impulses at an energy flux density of 0.41 mJ/mm² and a frequency of 4 Hz^{6,12}.



Postoperative Treatment

Postoperatively, a below-elbow thermoplastic splint incorporating the thumb up to the interphalangeal joint was applied for eight weeks. After removal of the splint, active and passive wrist exercises were initiated in conjunction with hand therapy.

Follow-up assessments included radiographs at the time of splint removal and CT scans at three and six months postoperatively. Upon radiographic confirmation of bony healing, strengthening exercises were initiated.

The CT scan performed six months after surgery (Figure 3a, 3b) demonstrated bony healing of the scaphoid nonunion. Due to an impingement between the scaphoid plate and the palmar rim of the distal radius, removal of the plate was planned.



Figure 3: CT images at 6 months follow-up



Conclusion

Recalcitrant scaphoid nonunions remain a substantial challenge for hand surgeons due to their limited biological healing potential. Consequently, some authors advocate salvage procedures rather than reconstructive strategies in these cases. However, angular-stable plate fixation has emerged as a reliable treatment option, demonstrating high union rates and favorable clinical outcomes¹⁻⁷. Compared with stabilization using only a single headless bone screw, plate fixation provides significantly greater rotational stability^{6,9}, and maintains space for cancellous bone graft insertion⁹. These aspects may contribute to the improved union rates reported in the literature¹⁻⁷. The availability of scaphoid plates in different sizes allows accommodation of the considerable anatomical variability of the scaphoid as well as the diversity of nonunion defects.

However, it should be considered that plate fixation is frequently associated with the need for secondary plate removal after nonunion healing¹⁻⁷, particularly when the plate needs to be positioned proximal to the so-called „equatorial-line“ to ensure sufficient stability². In such cases, impingement with the palmar rim of the distal radius may occur². Nevertheless, this should not be regarded as a complication but rather as an integral part of the treatment concept, and patients should be informed accordingly prior to surgery.¹⁰



References

1. Quadlbauer S, Pezzei C, Hintringer W, et al. Palmar Plate Fixation in Scaphoid Nonunions—Surgical Technique and Review of the Literature. *J Wrist Surg.* 2024;15(1):8-23. Published 2024 Nov 20. doi:10.1055/s-0044-1793941
2. Dodds SD, Williams JB, Seiter M, Chen C. Lessons learned from volar plate fixation of scaphoid fracture nonunions. *The Journal of hand surgery, European volume.* 2018/1// 2018;43(1):57–65. doi:10.1177/1753193417743636
3. Leixnering M, Pezzei C, Weninger P, et al. First Experiences With a New Adjustable Plate for Osteosynthesis of Scaphoid Nonunions. *J Trauma.* 2011/3// 2011;(1529-8809 (Electronic))
4. Putnam JG, DiGiovanni RM, Mitchell SM, Castañeda P, Edwards SG. Plate Fixation With Cancellous Graft for Scaphoid Nonunion With Avascular Necrosis. *Journal of Hand Surgery.* 2019;44(4):339.e1–339.e7. doi:10.1016/j.jhsa.2018.06.024
5. Putnam JG, Mitchell SM, DiGiovanni RM, Stockwell EL, Edwards SG. Outcomes of Unstable Scaphoid Nonunion With Segmental Defect Treated With Plate Fixation and Autogenous Cancellous Graft. *Journal of Hand Surgery.* 2019;44(2):160.e1–160.e7. doi:10.1016/j.jhsa.2018.05.023
6. Quadlbauer S, Pezzei C, Beer T, et al. Treatment of scaphoid waist nonunion by one, two headless compression screws or plate with or without additional extracorporeal shockwave therapy. *Archives of orthopaedic and trauma surgery.* 2019/2// 2019;139(2):281–293. doi:10.1007/s00402-018-3087-6
7. Schormans P, Brink P, Poeze M, Hannemann P. Angular Stable Miniplate Fixation of Chronic Unstable Scaphoid Nonunion. *Journal of Wrist Surgery.* 2018;07(01):024–030. doi:10.1055/s-0037-1603202
8. Schmidle G, Ebner HL, Klauser AS, Fritz J, Arora R, Gabl M. Correlation of CT imaging and histology to guide bone graft selection in scaphoid non-union surgery. *Arch Orthop Trauma Surg.* Oct 2018;138(10):1395–1405. doi:10.1007/s00402-018-2983-0
9. Jurkowitsch J, Dall'Ara E, Quadlbauer S, et al. Rotational stability in screw-fixed scaphoid fractures compared to plate-fixed scaphoid fractures. *Archives of orthopaedic and trauma surgery.* 2016/11// 2016;136(11):1623–1628. doi:10.1007/s00402-016-2556-z
10. Schormans PMJ, Van der Heijden ARY, Poeze M, Ten Bosch JA, Hannemann PFW. The influence of plate removal on functional and patient reported outcomes after scaphoid nonunion surgery. *J Hand Surg Eur Vol.* Jan 15 2026;17531934251413137. doi:10.1177/17531934251413137
11. Quadlbauer S, Pezzei C, Jurkowitsch J, et al. [Palmar angular stable plate fixation of nonunions and comminuted fractures of the scaphoid]. *Oper Orthop Traumatol.* 2019;31(5):433-446. doi:10.1007/s00064-019-00623-0
12. Quadlbauer S, Pezzei C, Jurkowitsch J, et al. Double screw versus angular stable plate fixation of scaphoid waist nonunions in combination with intraoperative extracorporeal shockwave therapy (ESWT). *Arch Orthop Trauma Surg.* 2023;143(7):4565-4574. doi:10.1007/s00402-023-04806-0

Disclaimer: This case report presents the outcome of an individual patient case and it does not imply any guarantee or warranty in regard to treatment success. A surgeon must always rely on her or his own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. The professional must always comply with the individual product's Instructions For Use (www.medartis.com/documentation/instructions-for-use) as well as all laws and regulations. Medartis is not giving any medical advice. The devices may not be available in all countries due to registration and/or medical practices. All content – such as texts, video and pictures – was created by healthcare professional mentioned in the case report. For further questions, please contact your Medartis representative (www.medartis.com). This information contains products with CE and/or UKCA marking.
For US only: Federal law restricts this device to sale by or on the order of a physician.

© Medartis 2025. Everything herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Medartis or its affiliates unless otherwise indicated. It is forbidden to redistribute, duplicate or disclose, anything herein, in whole or in part, without the prior written consent of Medartis.