


medartis

PRECISION IN FIXATION

PRODUCT INFORMATION

# TTS – Titanium Trauma Splint



# MODUS

# TTS – Titanium Trauma Splint



## Easy adaptation, fixation and removal

- Physiologic
- Small dimension
- Flexibility
- High comfort
- Cleansability

## Features

### For the dentist

- User-friendly handling
- Easy adaptation and fixation
- Manual shaping of the splint
- Easy to cut to size
- Time reduction in simple splint application and removal
- Optimized splint design

## Indications

The TTS Tooth Splint is used for:

- Tooth stabilization after repositioning or replantation (acute trauma cases)
- Tooth stabilization after intentional (planned) replantation

### For the patient

- Improved comfort: minimal impairment of speech, eating and oral hygiene in comparison with conventional splints
- Preservation of physiological mobility of the injured tooth as well as the neighboring teeth

## Application



Repositioning



Cleansing



Etching



Drying



Bonding agent



Composite



# Clinical Cases

## Case 1 (Michael Bornstein)



This 15-year-old patient was admitted following a road accident. Avulsion of 11 (21, 11) showed an enamel-dentine fracture without involvement of the pulp. 1 (21) was replanted and stabilized for 10 days by means of the TTS.

## Case 2 (Yves Germanier)



Complex tooth trauma in a 13-year-old girl after a fall while climbing. Avulsion of 1 (21), root fracture of 1 in the cervical third of the root (crown fragment could not be found any more) and 2 (12) with a fracture of the crown without involvement of the pulp. 1 (21) was replanted and fixated for two weeks using the TTS.

Developed in cooperation with: Department of Oral Surgery and Stomatology, School of Dental Medicine University of Berne, Switzerland, T. von Arx and colleagues. Cases and pictures received from the above.

## Publications

von Arx, T., Filippi, A., Buser, D. (2001a)  
**Splinting of traumatized teeth with a new device: TTS (Titanium Trauma Splint)**  
Dent Traumatol, 17(4), pp. 180-4.

von Arx, T., Filippi, A., Lussi, A. (2001b)  
**Comparison of a new dental trauma splint device (TTS) with three commonly used splinting techniques**  
Dent Traumatol, 17(6), pp. 266-74.

Filippi, A., von Arx, T., Lussi, A. (2002)  
**Comfort and discomfort of dental trauma splints - a comparison of a new device (TTS) with three commonly used splinting techniques**  
Dent Traumatol, 18(5), pp. 275-80.

Ingimarsson, S., von Arx, T. (2002)  
**Neue Schienungstechnik in der Zahntraumatologie  
Nouvelle technique de contention en traumatologie dentaire**  
[A new splint technique in dental traumatology]  
Schweiz Monatsschr Zahnmed, 112(12), pp. 1263-73. Article in German and French.

von Arx, T. (2005)  
**Splinting of traumatized teeth with focus on adhesive techniques**  
J Calif Dent Assoc, 33(5), pp. 409-14.

Chappuis, V., von Arx, T. (2005)  
**Replantation of 45 avulsed permanent teeth: a 1-year follow-up study**  
Dent Traumatol, 21(5), pp. 289-96.

Adatia, A., Kenny, D. J. (2006)  
**Titanium Trauma Splint: An Alternative Splinting Product**  
JCDA, 72(No. 8), pp. 721-23.

Ferrazzini Pozzi, E. C., von Arx, T. (2008)  
**Pulp and periodontal healing of laterally luxated permanent teeth: results after 4 years**  
Dent Traumatol, 24(6), pp. 658-62.

R\_TTS-01000001\_v0 / © 2021-03, Medartis AG, Switzerland. All technical data subject to alteration.

#### **MANUFACTURER & HEADQUARTERS**

Medartis AG | Hochbergerstrasse 60E | 4057 Basel/Switzerland  
P +41 61 633 34 34 | F +41 61 633 34 00 | [www.medartis.com](http://www.medartis.com)

#### **SUBSIDIARIES**

Australia | Austria | Brazil | China | France | Germany | Japan | Mexico | New Zealand | Poland | Spain | UK | USA

For detailed information regarding our subsidiaries and distributors, please visit [www.medartis.com](http://www.medartis.com)



Disclaimer: This information is intended to demonstrate the Medartis portfolio of medical devices. A surgeon must always rely on her or his own professional clinical judgement when deciding whether to use a particular product when treating a particular patient. Medartis is not giving any medical advice. The devices may not be available in all countries due to registration and/or medical practices. For further questions, please contact your Medartis representative ([www.medartis.com](http://www.medartis.com)). This information contains CE-marked products.  
For US only: Federal law restricts this device to sale by or on the order of a physician.