

medartis

PRECISION IN FIXATION

Hand Product Overview



APTUS

Contents

- 3 Introduction Medartis
- 4 Comprehensive Solutions for the Hand
- 6 APTUS Technologies
- 8 Plate Designs
 - 10 1.2/1.5 Fixation Plates
 - 11 1.2/1.5 TriLock Plates
 - 12 2.0/2.3 Fixation Plates
 - 13 2.0/2.3 TriLock Plates
 - 14 2.0/2.3 MC Compression Plates
 - 15 2.0/2.3 TriLock Arthrodesis Plates
 - 16 Special Plates
 - Hook
 - Scaphoid
 - Rotation (1.5 and 2.0)
 - Condylar (1.2/1.5 and 2.0/2.3)
- 18 Screws and Washers
- 19 Cannulated Compression Screws: CCS and headed CCS 1.7, 2.2, 3.0
- 20 Instruments
- 22 Bibliography
- 23 Loan Service and Contact Addresses

For further information regarding the APTUS product line visit
www.medartis.com.

Introduction Medartis

Medartis, headquartered in Basel, Switzerland, specializes in technical high-precision implant systems for surgical fixation of bone fractures and osteotomies.

Medartis develops, manufactures and sells titanium screws and plates, surgical instruments and system solutions for fracture fixation in the facial skull and the extremities.

Our motto is «Precision in fixation». Since the company's founding in 1997, we place the highest priority on maintaining stringent quality standards, continuous further development and innovation as well as comprehensive service provision.


Medartis is represented worldwide through its subsidiaries and a broad distributor network.




Comprehensive Solutions for the Hand

APTUS Hand

A broad portfolio of hand plates and screws to support the treatment of all fracture types and corrections.





APTUS Cannulated Compression Screws



Comprehensive portfolio of cannulated screws and headed cannulated screws available with short thread, long thread or fully threaded in different diameters and lengths.

APTUS Technologies

TriLock locking technology



All APTUS systems are based on the multidirectional and angular stable TriLock locking technology.

- Patented TriLock locking technology
 - multidirectional locking of the screw in the plate
 - spherical three-point wedge-locking
 - friction locking through radial bracing of the screw head in the plate without additional tensioning components
- Screws can pivot freely by $\pm 15^\circ$ in all directions for optimal positioning
- Fine-tuning capabilities of fracture fragments
- TriLock screws can be relocked in the same plate hole at individual angles up to three times
- Minimal screw head protrusion thanks to internal locking contour
- No cold welding between plate and screws




Patented TriLock locking technology – multidirectional locking of the screw in the plate

Minimal screw head protrusion



Screw technology

- Patented HexaDrive screw head design
 - HexaDrive interface with self-holding properties between screw and screwdriver
 - increased torque transmission
 - simplified screw pick-up due to patented self-holding technology
- Soft tissue protection due to smooth screw head design
- Atraumatic screw tip offers soft tissue protection when inserting screws bicortically
- Increased torsional, bending and shear stability due to conical core
- Precision-cut thread profile for sharpness and self-tapping properties
- Adapted pitches of the screw thread in cortical screws depending on screw length
- Double-threaded TriLock screws reduce screw insertion time



SpeedTip thread design

- Functionally unique cutting with immediate bite¹³
- Immediate cutting of the bone with only slight axial pressure
- The triangular tip design permits simultaneous drilling, tapping and compression of the bone tissue during insertion for increased pull-out stability^{14, 15}
- Reduced insertion torque thanks to the polygonal tip and tapered shaft




Plate Designs



General features and benefits

- Anatomically precontoured and fracture-specific implant geometries for less or no intraoperative bending
- Rounded edges and a smooth surface for soft tissue protection
- Offset screw holes in numerous plates avoid screw collisions



- Color-coded implants for easy identification in the operating room:
 - gold = Fixation/compression plates and cortical screws (fixation)
 - blue = TriLock plates and TriLock screws (locking)
- Two screw sizes can be used for each plate thickness:
 - 1.2/1.5 screws for 1.2/1.5 plates with a plate thickness of 0.6 or 0.8 mm
 - 2.0/2.3 screws for 2.0/2.3 plates with a plate thickness of 1.0 or 1.3 mm

1.2 / 1.5 Fixation Plates

Typical clinical findings


- Fractures in phalanges
- Unstable, non-comminuted fractures
- Simple intraarticular and extraarticular fractures
- For PIP arthrodesis

Features and benefits

- Several design options: straight, L-/T-/Y-shape, and grid structure
- Low plate profile with 0.6 mm thickness
- L-, Y- and T-plates
 - reduction of periarticular and intraarticular fractures
- Broader bars in bridging zones
- Grid plates
 - increased torsional and bending stability



Clinical Examples



Preoperative X-rays




Intraoperative image



Postoperative X-rays


Clinical example published with the kind permission of: C. Ranft, Kiel, Germany



Preoperative X-rays



Intraoperative image



Postoperative X-rays

Clinical example published with the kind permission of the author


1.2 / 1.5 TriLock Plates

Typical clinical findings


- Fractures in phalanges
- Comminuted fractures
- Simple and complex intraarticular and extraarticular fractures
- Fractures in osteoporotic bone
- For PIP arthrodesis

Features and benefits

- Internal fixator principle for early mobilization
- Several design options: straight, T-shape, and grid structure
- Low plate profile with 0.8 mm thickness
- Double bars between screw holes for increased torsional stability
- Double row T-plate:
 - more screw options for support in the subchondral area
 - reduction of articular fracture fragments and fractures close to the joint
- Grid plates:
 - increased torsional and bending stability
 - for comminuted fractures in the metaphyseal and diaphyseal region
 - straight hole arrangement on one end of the plate to position the plate as close to the joint as possible



Clinical Example



Preoperative X-rays



Intraoperative image



Postoperative X-rays

Clinical example published with the kind permission of: D. Schäfer, Basel, Switzerland

2.0 / 2.3 Fixation Plates

Typical clinical findings

- Fractures in metacarpals and proximal phalanges
- Unstable, non-comminuted fractures
- Simple intraarticular and extraarticular fractures

Features and benefits

- Several design options: straight, L-/T-/Y-shape, and grid structure
- Low plate profile with 1.0 mm thickness
- L-, Y- and T-plates
 - reduction of periarticular and intraarticular fractures
- Broader bars in bridging zones
- Grid plates
 - increased torsional and bending stability





2.0 / 2.3 TriLock Plates

Typical clinical findings


- Fractures in metacarpals and proximal phalanges
- Comminuted fractures
- Simple and complex intraarticular and extraarticular fractures
- Fractures in osteoporotic bone

Features and benefits

- Internal fixator principle for early mobilization
- Several design options: straight, L-/T-/Y-shape, and grid structure
- Low plate profile with either 1.0 or 1.3 mm thickness
 - choice of construct stability
- L-, Y- and T-plates:
 - reduction of articular fracture fragments and fractures close to the joint
 - broader bars in bridging zones
- Grid plates:
 - increased torsional and bending stability
 - for comminuted fractures in the metaphyseal and diaphyseal region
 - straight and offset hole arrangements



Clinical Example



Preoperative X-ray

Intraoperative images

Postoperative X-ray

Clinical example published with the kind permission of: C. Ranft, Kiel, Germany


2.0 / 2.3 MC Compression Plates

Typical clinical findings

- Fractures in metacarpals only
- Extraarticular fractures
- Not for comminuted fractures
- Transverse fractures (compression plating)
- Oblique and spiral fractures (neutral plating)

Features and benefits

- Several design options: straight, L-/Y-shape
- Low plate profile with 1.3 mm thickness
- Up to 2 mm compression possible
- Y- and T-plates
 - reduction of periarticular and intraarticular fractures
- Broader bars in bridging zones



Clinical Example



Preoperative X-rays

Intraoperative image


Postoperative X-rays

Clinical example published with the kind permission of the author

2.0 / 2.3 TriLock Arthrodesis Plates

Typical clinical findings

- 4CF plates
 - Osteoarthritis between radius, scaphoid and potentially midcarpal joint
- STT plate
 - Osteoarthritis between scaphoid, trapezium and trapezoid




Features and benefits

- Internal fixator principle for early mobilization
- Options for compression and stable fixation of the carpal bones
 - compression using cortical screws
 - angular stable fixation using TriLock screws
- Concave shape of plates and corresponding reamers
 - minimal bone removal
 - simple handling without intraoperative tilting of the plates
 - at least two screws can be inserted in each carpal bone
- 4CF-plates
 - two sizes to meet individual anatomical requirements
 - straight and offset hole arrangements




Clinical Examples



Intraoperative X-rays



Intraoperative image




Postoperative X-rays

Clinical example published with the kind permission of: Arnold-Peter Weiss, Providence, USA



Preoperative X-rays



Intraoperative image



Postoperative X-rays

Clinical case published with the kind permission of: Arnold-Peter Weiss, Providence, USA

Special Plates

Hook plate

Typical clinical findings

- Bony avulsions (e.g. mallet fractures)

Features and benefits

- Tap on plate for supported fragment reduction with holding instrument
- Compression hole



Scaphoid plate

Typical clinical findings

- Non-unions and multifragmentary fractures of the scaphoid

Features and benefits

- High degree of stability due to grid structure and up to three TriLock screws on each side of the non-union
- Anatomically preshaped to correct humpback deformity
- Two middle bars to keep the bone graft in place




Rotation correction plates

Typical clinical findings

- For rotational corrections in phalanges and metacarpals

Features and benefits

- Transversal oblong hole allows up to $\pm 25^\circ$ of rotational correction and is close to the joint to perform the osteotomy near the metaphyseal area
- Separate screw holes ("frog design") simplify contouring in the periarticular area




Condylar plates

Typical clinical findings


- For lateral reduction of simple intraarticular fractures

Features and benefits


- Lateral reduction of simple intraarticular fractures
- Compression holes




Clinical Examples



Preoperative X-ray




Intraoperative image




Postoperative X-ray


Clinical example published with the kind permission of: M. Aerni, Münsterlingen, Switzerland



Preoperative X-rays




Intraoperative images



Postoperative X-rays


Clinical example published with the kind permission of the author



Preoperative X-ray



Intraoperative image



Postoperative X-ray

Clinical example published with the kind permission of the author

Screws and Washers

Screws

– Choice of different screw diameter for each plate type:

	Cortical Screws				TriLock Screws	
Diameter (mm)	1.2	1.5	2.0	2.3	1.5	2.0
Plate types	1.2/1.5 Fixation Plates 1.2/1.5 TriLock Plates		2.0/2.3 Fixation Plates 2.0/2.3 TriLock Plates 2.0/2.3 MC Compression Plates 2.0/2.3 TriLock Arthrodesis Plates		1.2/1.5 TriLock Plates	2.0/2.3 TriLock Plates 2.0/2.3 TriLock Arthrodesis Plates

- Intraoperative adjustability of screw angle and fine-tuning capabilities
- Cortical screw for "screw-only" fracture fixation
- Emergency (cortical) screws available


Washers

- Biconcave design
- For force distribution over a larger bone surface around the screw head.

Cannulated Compression Screws: CCS and headed CCS 1.7, 2.2, 3.0

Typical clinical findings

- Head / base / oblique fractures in phalanges
- DIP joint arthritis
- Transverse fracture in metacarpals
- Bennett fracture
- Scaphoid fracture
- Ulnar styloid fracture
- Radial styloid fracture




Features and benefits

- Three different diameters: 1.7 mm, 2.2 mm and 3.0 mm
- Two different thread lengths
- Headed and headless screws
- Covering numerous areas of use in the hand
- Patented SpeedTip thread design
 - functionally unique cutting with immediate bite¹³
 - immediate cutting of the bone with only slight axial pressure
- The triangular tip design permits simultaneous drilling, tapping and compression of the bone tissue during insertion for increased pull-out stability^{14, 15}
- Reduced insertion torque thanks to the polygonal tip and tapered shaft


Clinical Examples



Preoperative X-ray



Intraoperative images





Postoperative X-rays


Clinical example published with the kind permission of: N. Schelhorn, R. Fricker, Switzerland


Instruments

- Simple and easy to use
- Intuitive application
- Consistent and distinct color coding:

APTUS 1.2 = red 

APTUS 1.5 = green 

APTUS 2.0 = blue 

APTUS 2.3 = brown 

Ergonomically designed and compact instrument kit

One plate holding and positioning instrument for each of the systems 1.2/1.5 and 2.0/2.3



One screwdriver/handle with screwdriver blade for each of the systems 1.2/1.5 and 2.0/2.3



Twist drills for core and gliding hole for all screw sizes 1.2, 1.5, 2.0 and 2.3



One drill guide for each of the systems 1.2/1.5 and 2.0/2.3



Specific drill guides for lag screw technique



Countersinks for 1.2/1.5 and 2.0/2.3 screws



Depth gauge for one-handed use



One pair of plate bending pliers for all plate sizes



One pair of plate cutting pliers for all plate sizes



Bibliography

1. Park, J.H., et al., Dorsal transosseous reduction and locking plate fixation for articular depressed middle phalangeal base fracture. *Arch Orthop Trauma Surg*, 2019. 139(1): p. 141-145.
2. Doering, T.A., A.S. Greenberg, and D.V. Tuckman, Dorsal Plating for Intra-articular Middle Phalangeal Base Fractures With Volar Instability. *Hand (N Y)*, 2019. 14(5): p. 620-625.
3. Vester, H., et al., Patient satisfaction after hook plate treatment of bony avulsion fracture of the distal phalanges. *Eur J Med Res*, 2018. 23(1): p. 35.
4. Spies, C.K., et al., [Open reduction and screw/plate osteosynthesis of metacarpal fractures]. *Oper Orthop Traumatol*, 2019. 31(5): p. 422-432.
5. Eng, K., et al., Volar Scaphoid Plating for Nonunion: A Multicenter Case Series Study. *J Wrist Surg*, 2020. 9(3): p. 225-229.
6. Schormans, P.M.J., et al., Mid-term outcome of volar plate fixation for scaphoid nonunion. *Bone Joint J*, 2020. 102-b(12): p. 1697-1702.
7. Putnam, J.G., et al., Outcomes of Unstable Scaphoid Nonunion With Segmental Defect Treated With Plate Fixation and Autogenous Cancellous Graft. *J Hand Surg Am*, 2019b. 44(2): p. 160.e1-160.e7.
8. Athlani, L., et al., Four-corner arthrodesis with a dorsal locking plate: 4-9-year follow-up. *J Hand Surg Eur Vol*, 2020. 45(7): p. 673-678.
9. Chaudhry, T., et al., Four corner fusion using a multidirectional angular stable locking plate. *World J Orthop*, 2016. 7(8): p. 501-6.
10. Giesen, T., et al., Intramedullary headless screw fixation for fractures of the proximal and middle phalanges in the digits of the hand: a review of 31 consecutive fractures. *J Hand Surg Eur Vol*, 2016. 41(7): p. 688-94.
11. Poggetti, A., et al., Outcomes of 173 metacarpal and phalangeal fractures treated by intramedullary headless screw fixation with a 4-year follow-up. *J Hand Surg Eur Vol*, 2021. 46(5): p. 466-470.
12. Ek, E.T., et al., Clinical Outcomes of Double-Screw Fixation with Autologous Bone Grafting for Unstable Scaphoid Delayed or Nonunions with Cavitory Bone Loss. *J Wrist Surg*, 2020(EFirst).
13. Spiegel, A.; Pochlatko, N.; Zeuner, H.; Lang, A.; Biomechanical Tests of Different Cannulated Compression Screws (on file; Medartis AG, Switzerland)
14. Heidemann, W.; Terheyden, H.; Gerlach, K. L.: Analysis of the osseous/metal interface of drill free screws and self-tapping screws (*Journal of Cranio-Maxillofacial Surgery*, 2001, 29, 69 – 74)
15. Heidemann, W.; Terheyden, H.; Gerlach, K. L.: In-vivo-Untersuchungen zum Schrauben-Knochen-Kontakt von Drill-Free-Schrauben und herkömmlichen selbstschneidenden Schrauben (*Mund Kiefer GesichtsChir* 5 2001: 17 – 21)

Medartis Loan Service and Contact Addresses

All APTUS systems are also available as a loan set:

24 hrs service (Monday–Friday): order today for delivery on the next working day*

Contact addresses

Contact us for further information regarding the APTUS product line:

Medartis Headquarters Switzerland	Phone: +41 (0)61 633 34 34 Fax: +41 (0)61 633 34 00 order@medartis.com
Australia	Phone: 1300 858 853/+61 (0)7 3326 8700 Fax: +61 (0)7 3862 2665 bookings@medartis.com
Austria	Phone: +43 (0) 5577 62 776 Fax: +43 (0) 5577 62 776 20 orders_at@medartis.com
Brazil	Phone: + 55 11 3624-7844 atendimento.br@medartis.com
France	Phone: +33 (0) 4 74 99 94 14 Fax: +33 (0) 4 74 99 00 19 commandes-fr@medartis.com
Germany	Phone: +49 (0) 7665 98 24 299 (loan service) Phone: +49 (0) 7665 98 24 0 Fax: +49 (0) 7665 98 24 10 orders_de@medartis.com
Japan	Phone: +81 3 4520 5048 Fax: +81 50 3737 5397 orders_jp@medartis.com
Mexico	Phone: (+52) 55 6388 7063 servicioclientes@medartis.com
New Zealand	Phone: 0800 548 001/+64 (9) 909 0416 Fax: +0800 548 002/+64 9 909 0419 bookings@medartis.com
Spain	Phone: +34 931446087 info.es@medartis.com
Poland	Phone: +48 (0) 71 359 56 18 Fax: +48 (0) 71 359 56 15 orders_pl@medartis.com
UK	Phone: +44 (0) 1924 47 66 99 Fax: +44 (0) 1924 47 20 00 orders_uk@medartis.com
USA	Phone: +1 574 376 2404 Fax: +1 574 966 1396 Toll Free: 877 406 BONE (2663) orders_us@medartis.com

* May vary by country



R_HAND-00001101_v0 / © 2022-11, Medartis AG, Switzerland.
All technical data subject to alteration.

MANUFACTURER & HEADQUARTERS

Medartis AG | Hochbergerstrasse 60E | 4057 Basel / Switzerland
P +41 61 633 34 34 | F +41 61 633 34 00 | www.medartis.com

SUBSIDIARIES

Australia | Austria | Brazil | France | Germany | Japan | Mexico |
New Zealand | Poland | Spain | UK | USA

For detailed information regarding our subsidiaries and distributors,
please visit www.medartis.com



Disclaimer: This information is intended to demonstrate the Medartis portfolio of medical devices. A surgeon must always rely on her or his own professional clinical judgement when deciding whether to use a particular product when treating a particular patient. Medartis is not giving any medical advice. The devices may not be available in all countries due to registration and/or medical practices. For further questions, please contact your Medartis representative (www.medartis.com). This information contains products with CE and/or UKCA marking. All pictures shown are for illustration purposes only and may not be an exact representation of the product.

For US only: Federal law restricts this device to sale by or on the order of a physician.