

Hand 1.2/1.5, 2.0/2.3



APTUS® Hand

Fixation Plates Hand 1.2/1.5



Scale 1 : 1

Order quantity: 1

Plate thickness (mm)




Fixation Plates Hand 2.0/2.3





MC (Metacarpal) Compression Plates Hand 2.0/2.3

4CF Plate 2.0/2.3



Locking Plates Hand 2.0



Instrument Case and Instruments


A-6020
Instrument Case with Silicon Inlay

A-6025
Instrument Tray

- A-2040**
Plate Bending Pliers with Vario Pin
- A-2060**
Plate and Screw Holding Forceps, angled
- A-2041**
Vario Plate Cutting Pliers for Hand Plates



- A-2610**
Screwdriver 2.0/2.3, HexaDrive 6
- A-2030**
Depth Gauge 1.2 - 2.3
- A-2650**
Plate Holding and Positioning Instrument 2.0/2.3



- A-2025**
Drill Guide 1.2 - 1.5
- A-2020**
Drill Guide 1.2 - 2.3
- A-2310**
Screwdriver 1.2/1.5, HexaDrive 4
- A-2350**
Plate Holding and Positioning Instrument 1.2/1.5

Screws, Washers, Twist Drills

Cortical Screws – Ø 1.2, HexaDrive 4



Art. No.	Length	Pieces per Pack
A-5100.04	4 mm	5
A-5100.05	5 mm	5
A-5100.06	6 mm	5
A-5100.07	7 mm	5
A-5100.08	8 mm	5
A-5100.09	9 mm	5
A-5100.10	10 mm	5
A-5100.11	11 mm	5
A-5100.12	12 mm	5
A-5100.13	13 mm	5
A-5100.14	14 mm	5
A-5100.16	16 mm	5
A-5100.18	18 mm	5
A-5100.20	20 mm	5

For a package of 1 please add '11' at the end of the art. no., for example A-5100.2011


Cortical Screws – Ø 1.5, HexaDrive 4



Art. No.	Length	Pieces per Pack
A-5200.04	4 mm	5
A-5200.05	5 mm	5
A-5200.06	6 mm	5
A-5200.07	7 mm	5
A-5200.08	8 mm	5
A-5200.09	9 mm	5
A-5200.10	10 mm	5
A-5200.11	11 mm	5
A-5200.12	12 mm	5
A-5200.13	13 mm	5
A-5200.14	14 mm	5
A-5200.15	15 mm	5
A-5200.16	16 mm	5
A-5200.17	17 mm	5
A-5200.18	18 mm	5
A-5200.19	19 mm	5
A-5200.20	20 mm	5
A-5200.21	21 mm	5
A-5200.22	22 mm	5
A-5200.23	23 mm	5
A-5200.24	24 mm	5

For a package of 1 please add '11' at the end of the art. no., for example A-5200.2011


Emergency Cortical Screws – Ø 1.8, HexaDrive 4



Art. No.	Length	Pieces per Pack
A-5300.06	6 mm	5
A-5300.10	10 mm	5

For a package of 1 please add '11' at the end of the art. no., for example A-5300.1011

Biconcave Washers – Ø 1.2/1.5, 2.0/2.3



Art. No.	System	Pieces per Pack
A-4300.70	1.2/1.5	5
A-4600.70	2.0/2.3	5

1 : 1

TriLock Locking Screws – Ø 2.0, HexaDrive 6



Art. No.	Length	Pieces per Pack
A-5450.06	6 mm	5
A-5450.07	7 mm	5
A-5450.08	8 mm	5
A-5450.09	9 mm	5
A-5450.10	10 mm	5
A-5450.11	11 mm	5
A-5450.12	12 mm	5
A-5450.13	13 mm	5
A-5450.14	14 mm	5
A-5450.16	16 mm	5
A-5450.18	18 mm	5
A-5450.20	20 mm	5

For a package of 1 please add '11' at the end of the art. no., for example A-5450.2011

Cortical Screws – Ø 2.0, HexaDrive 6



Art. No.	Length	Pieces per Pack
A-5400.04	4 mm	5
A-5400.05	5 mm	5
A-5400.06	6 mm	5
A-5400.07	7 mm	5
A-5400.08	8 mm	5
A-5400.09	9 mm	5
A-5400.10	10 mm	5
A-5400.11	11 mm	5
A-5400.12	12 mm	5
A-5400.13	13 mm	5
A-5400.14	14 mm	5
A-5400.15	15 mm	5
A-5400.16	16 mm	5
A-5400.17	17 mm	5
A-5400.18	18 mm	5
A-5400.19	19 mm	5
A-5400.20	20 mm	5
A-5400.21	21 mm	5
A-5400.22	22 mm	5
A-5400.23	23 mm	5
A-5400.24	24 mm	5

For a package of 1 please add '11' at the end of the art. no., for example A-5400.1011

Cortical Screws – Ø 2.3, HexaDrive 6



Art. No.	Length	Pieces per Pack
A-5500.05	5 mm	5
A-5500.06	6 mm	5
A-5500.07	7 mm	5
A-5500.08	8 mm	5
A-5500.09	9 mm	5
A-5500.10	10 mm	5
A-5500.11	11 mm	5
A-5500.12	12 mm	5
A-5500.13	13 mm	5
A-5500.14	14 mm	5
A-5500.15	15 mm	5
A-5500.16	16 mm	5
A-5500.17	17 mm	5
A-5500.18	18 mm	5
A-5500.19	19 mm	5
A-5500.20	20 mm	5
A-5500.21	21 mm	5
A-5500.22	22 mm	5
A-5500.23	23 mm	5
A-5500.24	24 mm	5
A-5500.26	26 mm	5
A-5500.28	28 mm	5
A-5500.30	30 mm	5
A-5500.32	32 mm	5
A-5500.34	34 mm	5

For a package of 1 please add '11' at the end of the art. no., for example A-5500.1011

Emergency Cortical Screws – Ø 2.5, HexaDrive 6



Art. No.	Length	Pieces per Pack
A-5600.06	6 mm	5
A-5600.10	10 mm	5

For a package of 1 please add '11' at the end of the art. no., for example A-5600.1011

Twist Drills for 1.2 mm Screws



Art. No.	Ø	Stop	Length	Drill Shaft End	Pieces per Pack
A-3110	1.0 mm	20 mm	82 mm	Dental	1
A-3120	1.0 mm	20 mm	82 mm	Stryker J-Latch	1
A-3130	1.0 mm	20 mm	76 mm	AO Quick Coupling	1
A-3112	1.0 mm	20 mm	52 mm	Dental	1

Lag Screw Twist Drills for 1.2 mm Screws



Art. No.	Ø	Stop	Length	Drill Shaft End	Pieces per Pack
A-3111	1.2 mm	10 mm	72 mm	Dental	1
A-3121	1.2 mm	10 mm	72 mm	Stryker J-Latch	1
A-3131	1.2 mm	10 mm	66 mm	AO Quick Coupling	1
A-3113	1.2 mm	10 mm	42 mm	Dental	1

1 : 1

Twist Drills for 1.5 mm Screws



Art. No.	Ø	Stop	Length	Drill Shaft End	Pieces per Pack
A-3210	1.2 mm	25 mm	87 mm	Dental	1
A-3220	1.2 mm	25 mm	87 mm	Stryker J-Latch	1
A-3230	1.2 mm	25 mm	81 mm	AO Quick Coupling	1
A-3212	1.2 mm	25 mm	57 mm	Dental	1

Lag Screw Twist Drills for 1.5 mm Screws



Art. No.	Ø	Stop	Length	Drill Shaft End	Pieces per Pack
A-3211	1.6 mm	10 mm	72 mm	Dental	1
A-3221	1.6 mm	10 mm	72 mm	Stryker J-Latch	1
A-3231	1.6 mm	10 mm	66 mm	AO Quick Coupling	1
A-3213	1.6 mm	10 mm	42 mm	Dental	1

1 : 1

Twist Drills for 2.0 mm Screws



Art. No.	Ø	Stop	Length	Drill Shaft End	Pieces per Pack
A-3410	1.6 mm	25 mm	87 mm	Dental	1
A-3420	1.6 mm	25 mm	87 mm	Stryker J-Latch	1
A-3430	1.6 mm	25 mm	81 mm	AO Quick Coupling	1
A-3412	1.6 mm	25 mm	57 mm	Dental	1

Lag Screw Twist Drills for 2.0 mm Screws



Art. No.	Ø	Stop	Length	Drill Shaft End	Pieces per Pack
A-3411	2.1 mm	10 mm	72 mm	Dental	1
A-3421	2.1 mm	10 mm	72 mm	Stryker J-Latch	1
A-3431	2.1 mm	10 mm	66 mm	AO Quick Coupling	1
A-3413	2.1 mm	10 mm	42 mm	Dental	1

1 : 1

Twist Drills for 2.3 mm Screws



Art. No.	Ø	Stop	Length	Drill Shaft End	Pieces per Pack
A-3510	1.9 mm	35 mm	97 mm	Dental	1
A-3520	1.9 mm	35 mm	97 mm	Stryker J-Latch	1
A-3530	1.9 mm	35 mm	91 mm	AO Quick Coupling	1
A-3512	1.9 mm	25 mm	57 mm	Dental	1

Lag Screw Twist Drills for 2.3 mm Screws



Art. No.	Ø	Stop	Length	Drill Shaft End	Pieces per Pack
A-3511	2.35 mm	10 mm	72 mm	Dental	1
A-3521	2.35 mm	10 mm	72 mm	Stryker J-Latch	1
A-3531	2.35 mm	10 mm	66 mm	AO Quick Coupling	1
A-3513	2.35 mm	10 mm	42 mm	Dental	1

1 : 1

Lid for Instrument Case, Countersinks and Reamer

A-6022

Lid Hand 1.2 – 2.3 for Instrument Case A-6020



Countersinks for Lag Screws



Art. No.	System	Ø	Length	Countersink Shaft End	Pieces per Pack
A-3310	1.2/1.5	2.7 mm	37 mm	Dental	1
A-3610	2.0/2.3	3.5 mm	37 mm	Dental	1

1 : 1

Reamer Four Corner Fusion Plate



Art. No.	System	Ø	Length	Reamer Shaft End	Pieces per Pack
A-3630	2.0/2.3	17 mm	87 mm	AO Quick Coupling	1

1 : 2