

medartis®

PRECISION IN FIXATION

PRODUCT INFORMATION

# OptoFix 2.0 Transbuccal Micro Endoscope



MODUS®  
Mandible

# OptoFix 2.0

## Transbuccal Micro Endoscope

with Integrated Working Channel for Endoscopically Assisted Internal Fixation of Mandibular Fractures

### Product Features

- Simplified handling in difficult to access fractures
- Direct transbuccal endoscopic vision
- Integrated fibre optical illumination
- Direct drilling and screw insertion through working channel
- Integrated irrigation duct
- Instrumentation adapted for MODUS TriLock 2.0/2.3/2.5, MODUS 2.0 and MODUS Trauma 2.0



### Indications

- Endoscopic transbuccal osteosynthesis of fractures involving the condylar neck, condylar base and ramus region of the mandible
- Supports the intraoral, retromandibular and accordingly submandibular/angular endoscopically assisted approach

# Ordering Information

## Micro Endoscope OptoFix



## Protection Sleeve OptoFix



## Trocar OptoFix



## Adapter of M-2830



Art. No.	∅	Description	Length	Piece/Pkg
M-2830	4.7 mm		102 mm*	1
M-2840			34 mm	1
M-2850		for M-2840	41 mm	1
M-2830.1	11 mm	for Storz Light Guide	22 mm	1
M-2830.2	9 mm	for Wolf Light Guide	7 mm	1

\* Length of handle/working channel

### Soft Tissue Retraction Clamp



Art. No.	Description	Length	Piece/Pkg
M-2860	for M-2840	200 mm	1

### Screwdriver Blades 2.0



M-2823



M-2813

Art. No.	Interface	Description	Length	Piece/Pkg
M-2823	HD6	self-holding	130 mm	1
M-2813		self-holding	130 mm	1

### Twist Drills for Screws Ø 2.0 mm and Cancellous Screws Ø 2.3 mm

#### Ø 1.5 mm



M-3443



M-3433

Art. No.	Ø Twist Drill	Stop	Length	Drill Shaft End	Piece/Pkg
M-3443	1.5 mm	10 mm	141 mm	Dental	1
M-3433	1.5 mm	10 mm	141 mm	Stryker J-Latch	1

### Screwdriver Handle 2.0



1:1

Art. No.	Piece/Pkg
M-2503	1

### Instrument Case



M-6320

### Lid for Instrument Case



M-6330

Art. No.	Description	Size	Piece/Pkg
M-6320	steel	231 x 244 x 54 mm	1
M-6330	steel	228 x 243 x 10 mm	1

### Instrument Trays for M-6320



M-6350



M-6360

Art. No.	Description	Piece/Pkg
M-6350	for micro endoscope	1
M-6360	for instruments	1

### Coding Label



Art. No.	Description	Size	Piece/Pkg
M-6535	MODUS® OptoFix 2.0	20 x 60 mm	1

MANDIBLE-05000001\_v0 / © 01.2011, Medartis AG, Switzerland. All technical data subject to alteration.

#### **HEADQUARTERS**

Medartis AG | Hochbergerstrasse 60E | 4057 Basel/Switzerland

P +41 61 633 34 34 | F +41 61 633 34 00 | [www.medartis.com](http://www.medartis.com)

#### **SUBSIDIARIES**

Australia | Austria | France | Germany | Mexico | New Zealand | Poland | UK | USA

For detailed information regarding our subsidiaries and distributors, please visit [www.medartis.com](http://www.medartis.com)