

CMX PORTAL

| Date: 28. Jan. 2026 | Number of pages: 1

| Version: 2

WHITELISTING INSTRUCTIONS

The Medartis CMX Web Portal needs access to the following URLs in order to work properly:

URL	Ports	Schemes	Protocol
medartis.com/*	80, 443	http, https	TCP
cmx.medartis.com/*	80, 443	http, https, wss	TCP
medartiscustomerprod.b2clogin.com/*	80, 443	http, https	TCP
blbsigninui.blob.core.windows.net/*	80, 443	http, https	TCP
medartis-b2c-secure-ts.azurewebsites.net/*	80, 443	http, https	TCP